

RELAX AND RENEW

WITH HOT & COLD THERAPY GEL PACK

Aqua Bead Therapeutic Gel Pack Applies Heat Or Cold To Any Muscle.
 Use instructions are printed on the back of each piece. CPSIA, FDA and TRA Compliant.
 Reusable and Non-Toxic. Microwave & Freezer Safe.



Sideview



CPSIA, FDA and TRA Compliant

#GP100 Renew Hot & Cold Pack

Therapeutic Gel Pack applies heat or cold to sore muscles. Place in freezer for hot pack use and Microwave for hot pack use. Microwave and freezer safe. Non-toxic and reusable, safety tests available. Use instructions printed on the back of each piece. Extra large beads provide better performance and comfort.

Size: 4-1/2" Diameter

Imprint: 2-1/2" Diameter

150	250	500	1,000	2,500+	
1.60	1.45	1.35	1.25	1.19	(C)
	1.35	1.25	1.19		(C)

Available Colors:



Red



White



Blue



Lime Green



Use instructions printed on the back of each piece.

Screen Setup Charge: 55.00 (G). **Price Includes:** One color imprint, one location. One color imprint only. Standard non metallic. **Email Proof:** 7.50 (G), add 2 days to production time. **Repeat Setup Charge:** 35.00 (G). **Production Time:** 5 Days. Production time starts after all proof approvals, if required. (Rush service not available on this item). **Packaging:** Bulk. **Weight:** 34 lbs. (100pcs.)

Please specify Promo #E3291 to receive special pricing. Expires 8/31/2017. All Prices are in US Dollars

Note: After freezing, this item stays cold for up to 3 hours. If dropped, this item will not shatter or break.
 After microwaving 30 seconds, this item will remain warm up to 90 minutes. Condensation may potentially cause bubbling until cool.
 CAUTION: Product will be HOT after microwaving. Please refer to the Packaging Details table located on this page.
 0503

HAVE MORE TIME? TRY OVERSEAS PRODUCTION

Overseas Production:	5,000	10,000	25,000+
#GP100OP 8-10 Weeks	.85	.77	.69 (C)