

# QUICK-SHIP FULL COLOR SUBLIMATION FANDANA™ GAITER



\*Loose Fandana  
Shown



- Price Includes Full Color Imprint On Front Area Of Gaiter.
- All Purpose Unisex Accessory – Facemask, Headband, Beanie, Skullcap, Scrunchie, Neck Warmer, Muffler.
- One Size Fits All!

## #7134HT "THE FANDANA™ GAITER Multi-Functional Tubular Face, Head And Neck Wear

100	250	500	1,000	2,500+
3.00	2.90	2.80	2.70	2.59 (C)

**Price Includes:** A full color sublimation imprint on the front area of the Gaiter. **For additional location:** add .85 (G) running charge. **Setup Charge:** 60.00 (G). **Pre-Production Proofs:** 60.00 (G), 5 day lead time. **Repeat Setup Charge:** 27.50 (G). **Email Proof Charge:** 7.50 (G), add 2 day. **Material:** 100% Polyester Microfiber 130 GSM. **Production Time:** 5 Days. **Packaging:** Bulk. **Size:** 9-1/2" W x 18" H. **Imprint Area:** 6-1/2" W x 5-1/4" H. **Color:** White Only.

Specify Promo #E4198 to receive special pricing at time of order placement as it cannot be applied later. This offer may not be combined with any other promotion, discount, coupon and/or special column pricing (EQP & NQP do not apply sorry no exceptions). Expires 5/31/2021. All prices are in US Dollars.

This face cover is for source control use only. This face cover should not be used: (1) in any surgical setting or to provide liquid barrier protection; (2) in a clinical setting where the infection risk level through inhalation exposure is high; (3) for particulate filtration; or (4) in high risk aerosol generating procedures. The product has not been FDA cleared or approved, it has been authorized by the FDA under an EUA for use as source control by the general public during the COVID-19 pandemic. Source control means preventing the transmission of infection through a person's respiratory secretions which are produced when speaking, coughing, or sneezing. Face masks, including cloth face coverings, help with source control by covering the wearer's mouth and nose. This product is authorized only for the duration of the EUA during the COVID-19 outbreak, under section 564(b)(1) of the Act, 21 U.S.C. §360bbb-3(b)(1) unless the authorization is terminated or revoked sooner. If you experience any issues with this product, please contact us at [info@innovation-line.com](mailto:info@innovation-line.com). Important Information: Proofs provided are for logo placement and size approval only. The colors shown on e-proofs may not be accurately represented on your computer monitor (All metallic inks will NOT show up as metallic on screen) and PMS colors as printed in the Pantone Color Book are only a reference point. Therefore, e-proofs are not an exact color representation. Variations in print substrates/materials make it impossible to exactly duplicate PMS colors on Innovation Line products. Color variation cannot be considered a defect due to ink, lighting, and material variations. Innovation Line is not responsible for any color variation. In screen printing pad printing, Sublimation or transfer process exact color matches cannot be achieved. Regarding color matching, dye-sublimation printing uses a computer printer to generate a heat transfer to transfer the image onto fabric. Dye-sublimation printers use CMYK (Cyan Magenta Yellow Overcoating) colors, which differs from the more recognized CMYK colors in that the black is eliminated in favor of a clear overcoating. As a result, PMS inks cannot be matched exactly. They are only a reference point and an approximation 0219