


# FULL COLOR SUBLIMATION FANDITO GAITER

The new Fandito gaiter is half the size of our standard Fandana gaiter and designed primarily for use as a reusable facemask. An ultra-comfortable fit without annoying and uncomfortable ear loops. **Pricing includes up to a full color, full bleed, imprint on the outside (front and back) of the gaiter without any confusing additional color print charges. Pricing also includes freight delivery (in the continental US) to your customer's door in 12-15 days (door to door delivery), with no additional domestic freight cost or lag time. 130GSM thick polyester microfiber material. Individually polybagged.**

**DELIVERED TO YOUR  
CLIENT'S DOOR**  
(ANYWHERE IN THE CONTINENTAL US)  
**FREIGHT COST INCLUDED**



(Printed full color  
on the outside)



**#7144HTOP3 'THE FANDITO'**  
**MULTI-FUNCTIONAL FANDITO GAITER,  
TUBULAR HEAD & NECK WEAR**  
130GSM thick polyester microfiber material.

250	500	1,000	2,500+
2.40	2.10	1.90	1.75 (C)

**Size:** Approximately 9-3/4" W x 9-1/2" H.

**Imprint Area:** 9-3/4" W x 9-1/2" H, Full Bleed.

**Price Includes:** up to a full color, full bleed, imprint on the outside (front and back) of the gaiter. Pricing includes freight delivery (in the continental US) to your customer's door in 12-15 days (door to door delivery), with no additional domestic freight cost or lag time. Individually polybagged.

**Setup Charge:** 60.00 (G). **Pre-Production Proofs:** 60.00 (G), 5 day lead time. **Repeat Setup Charge:** 27.50 (G). **Email Proof Charge:** 7.50 (G), add 2 day. **Production Time:** 12-15 Days. **Color:** White Only.

Expires 5/31/2021.

This face cover is for source control use only. This face cover should not be used: (1) in any surgical setting or to provide liquid barrier protection; (2) in a clinical setting where the infection risk level through inhalation exposure is high; (3) for particulate filtration; or (4) in high risk aerosol generating procedures. The product has not been FDA cleared or approved, it has been authorized by the FDA under an EUA for use as source control by the general public during the COVID-19 pandemic. Source control means preventing the transmission of infection through a person's respiratory secretions which are produced when speaking, coughing, or sneezing. Face masks, including cloth face coverings, help with source control by covering the wearer's mouth and nose. This product is authorized only for the duration of the EUA during the COVID-19 outbreak, under section 564(b)(1) of the Act, 21 U.S.C. §360bbb-3(b)(1) unless the authorization is terminated or revoked sooner. If you experience any issues with this product, please contact us at [info@innovation-line.com](mailto:info@innovation-line.com).

Important Information: Proofs provided are for logo placement and size approval only. The colors shown on e-proofs may not be accurately represented on your computer monitor (All metallic inks will NOT show up as metallic on screen) and PMS colors as printed in the Pantone Color Book are only a reference point. Therefore, e-proofs are not an exact color representation. Variations in print substrates/materials make it impossible to exactly duplicate PMS colors on Innovation Line products. Color variation cannot be considered a defect due to ink, lighting, and material variations. Innovation Line is not responsible for any color variation. In screen printing pad printing, Sublimation or transfer process exact color matches cannot be achieved. Regarding color matching, dye-sublimation printing uses a computer printer to generate a heat transfer to transfer the image onto fabric. Dye-sublimation printers use CMYD (Cyan Magenta Yellow Overcoating) colors, which differs from the more recognized CMYK colors in that the black is eliminated in favor of a clear overcoating. As a result, PMS inks cannot be matched exactly. They are only a reference point and an approximation 0318